

Account Customization

Accessing Account Settings

The screenshot displays the Diagnocat 'Patients' interface. At the top, the Diagnocat logo is on the left, and navigation icons (home, globe, sun, and user profile) are on the right. The user profile icon is highlighted with a red box, indicating the path to access settings. Below the header, the 'Patients' section features a '+ Add new patient' button and a search bar labeled 'Search by Patient name or ID'. A grid of patient profiles is shown, each with a dental X-ray or photo, the patient's name, date of birth, last modified date, and the user's name (David Johnston). The user profile dropdown menu, also highlighted with a red box, shows the user's name 'David Johnston', a 'My Profile' link, the 'Current clinic' 'Dental Clinic' with a 'Clinic settings' link, and a 'Sign Out' button. A small robot icon is visible on the right side of the interface.

Diagnocat

David Johnston

Patients

+ Add new patient

All 14 Received Referrals 0 Sent Referrals 0

Search by Patient name or ID

Karen Morgan
DOB: July 8, 1971
Last modified: January 20, 2026
David Johnston

Abigail Foster
DOB: September 12, 1991
Last modified: January 20, 2026
KA

Nicole Schutz
DOB: December 18, 1980
Last modified: January 20, 2026
KA Katherine Amerson

Emma Smith
DOB: August 10, 1989
Last modified: January 20, 2026
KA

Joyce Morgan
DOB: August 14, 1972
Last modified: October 23, 2025
David Johnston

Karolin Lehmann
DOB: August 6, 2004
Last modified: October 23, 2025
David Johnston

Ethan Thomas
DOB: May 16, 1990
Last modified: October 23, 2025
David Johnston

Linda Bürger
DOB: June 9, 1976
Last modified: October 14, 2025
David Johnston

© Diagnocat 2026 Terms and conditions Data Processing Addendum The Privacy Policy About

David Johnston
My Profile

Current clinic
Dental Clinic
Clinic settings

Sign Out

1. Navigate to Settings: Click on your profile details in the top-right corner of the screen.
2. Choose Settings: Select either 'Clinic Settings' or 'Profile Settings' from the dropdown menu.

Profile Settings

In the 'My Profile' section, you can:

Update your personal information



Edit your first name, last name, and email address.

Add a profile photo

Upload a picture to personalize your account.

Upload your signature and stamp

On the right side, upload your signature and stamp in PNG format. These will automatically appear on your final patient reports, adding a professional touch.



Clinic Settings

Within 'Clinic Settings,' there are three tabs:

- 1

Clinic Information



 - Enter Clinic Details: Input your clinic's name and contact information.
 - Upload Branding Elements: Add your clinic's logo and official seal (PNG format) to ensure all AI-powered patient reports are still clinic branded and professional.
 - Choose dental notation format: Universal or ISO.
- 2

Team

 - Manage Team Access: Click 'Add Team Member' to invite new staff by entering their email, first name, last name, position, and assigning an access level.
 - Access Levels: Full Access: Complete platform access. Clinical Access: Access to clinical features. Limited Clinical Access: Restricted clinical features. Non-clinical Member: Access to administrative features only.
 - Edit or Remove Members: Use the icons on the right-hand side to update details or remove team members as needed.
- 3

Billing Information

 - Subscription Overview: Review your current subscription plan and usage statistics to monitor report usage at a glance.



Access levels

Access levels you can assign to team members

Access levels	Full access	Clinical access	Limited clinical access	Non-clinical Member
Manage company settings	+	-	-	-
Manage payments and subscriptions	+	-	-	-
View and manage team	+	+	-	-
View and manage patient's files	+	+	+	+
Share patient's profile	+	+	+	-
Edit patient's reports	+	+	+	-

By customizing these settings, Diagnocat becomes a tailored solution that fits seamlessly into your clinic's workflow.